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CAMSIZER® X2+

Shorter measurement time. Better optical resolution. More frames per second.

The CAMSIZER X2+ is one of the most powerful, extremely versatile particle size and shape analyzers with a wide measuring range that combines state-of-the-art camera technology with flexible dispersion options. Based on the principle of Dynamic Image Analysis (ISO 13322-2), the CAMSIZER X2+ provides precise particle size and shape information of powders, granules and suspensions in a measuring range from 0.9 μm to 8 mm.

The CAMSIZER X2+ produces a particle flow which is characterized by an optical system with high resolution. An ultrabright LED stroboscopic light sources and two high-resolution digital cameras achieve a frame rate of more than 420 images per second which are evaluated in real time by a powerful software. Thus, the CAMSIZER X2+ captures the images of hundreds of thousands to several millions of particles with highest accuracy within only 1 to 3 minutes.

The CAMSIZER X2+ provides a wide selection of particle information which allows for comprehensive and reliable characterization of the sample material. It is suitable for use in R&D as well as for routine tasks in quality control.

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

Many materials tend to agglomerate due to surface forces. The X-Jet module effectively disperses the particles in an air flow while passing through a Venturi nozzle. The dispersion pressure can be set as required for individual particles. For sensitive granules, for example, reduced pressure ensures non-destructive measurement. After passing the measurement field, the sample is automatically removed from the analyzer by a vacuum cleaner.

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
TRECHTERS EN GOTEN

To achieve optimum feeding conditions, even for adhesive powders, we offer different chutes made of either stainless steel or aluminium hard-coat. The hoppers can be made of aluminium or stainless steel and can hold up to 0.6 l of sample. The hopper is adjustable in height which promotes uniform flow of the sample.

The hoppers can be made of aluminium or stainless steel and can hold up to 0,6 l of sample. The hopper is adjustable in height which promotes uniform flow of the sample.

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

DEELTJES ANALYSER CAMSIZER X2+

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Due to the unique two-camera principle, all measurement ranges can be analyzed in one measurement run, without any hardware adjustments.

Measurement range specifications for the diameter of compact particles ¹⁾:


- | X-Jet (air pressure dispersion)
Standard: 0.9 μm to 2.5 mm with nozzle size 3.8 mm x 14 mm
Alternative: 0.9 μm to 6 mm with nozzle size 9 mm x 14 mm
Alternative: 5 μm to 8 mm with nozzle size 13.5 mm x 14 mm
- | X-Fall (gravity dispersion) ²⁾
Standard: 10 μm to 4 mm with nozzle size 7 mm x 14 mm
Alternative: 10 μm to 8 mm with nozzle size 14 mm x 14 mm
- | X-Flow (wet dispersion) ³⁾
Standard: 0.9 μm to 1 mm

1) Measurement ranges are sample dependent

2) The lower limit for X-Fall depends on the sample properties (agglomeration). For agglomerating particles the X-Jet or X-Flow module are recommended.

3) The upper size limit for X-Flow depends on the density of the sample material.

Measurement range specifications for fibers and elongated particles available on request.

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
TECHNISCHE GEGEVENS

Meetprincipe	Dynamische beeldanalyse (ISO 13322-2)
Meetbereik	0.9 µm to 8 mm 0.9 µm to 8 mm (air pressure dispersion) 10 µm to 8 mm (gravity dispersion) 0.9 µm to 1 mm (wet dispersion)
Type analyse	droge en natte analyse
Measurement time	~ 1 to 3 min (depends on desired measuring statistics, including analysis and data processing)
Aantal camera's	2 (Dual Camera Technology)
Monstervolume	< 20 mg - 500 g (afhankelijk van monstertype en meetwijze)
Camera performance	> 420 image/s, each with > 5.0 MPixel
Breedte van de meetzone	~ 350 mm ²
Digital Resolution	0.9 µm per pixel
Meetparameters	particle size (smallest diameter, length, mean diameter, etc.) particle shape (aspect ratio width to length, symmetry, sphericity, convexity etc., acc. to ISO 9276-6)
Afmetingen (B x H x D)	~ 850 x 580 x 570 mm
Gewicht (Meeteenheid)	~ 50 kg
Bedieningseenheid	Quad Core PC incl. Windows 11, monitor, keyboard and mouse, network card, PC interface cards for hardware communication, evaluation software

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