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MORE BELCONTROL FEATURES

- | Adsorption / desorption isotherms overlay and port comparison during measurement
- | Pressure, temperature and valve actuation trend data
- | System check for instrument diagnostics
- | E-mail notifications of measurement status and results
- | Interactive program in Japanese or English
- | Step-by-step instructions and extensive help functions

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
BELCONTROL: NEW OPERATION SOFTWARE

The versatility of BELSORP equipment is complemented by BELCONTROL, the intuitive, user-friendly operation software. It guides the user step-by-step through setup of analysis conditions, measurements, filling the liquid nitrogen or other bath, replacement of gas cylinders and degassing steps. It is designed to make the instrument accessible to inexperienced users.

For inexperienced users or measurements of unknown samples, BELCONTROL only requires basic sample information (name, mass, etc.), pretreatment conditions (if not performed externally) and measurement range.

Professional Mode permits detailed control of dosing settings, equilibrium criteria, leak-test options and other measurement configurations. Simple Mode reduces setup requirements; when a comparable prior isotherm is available, GDO can shorten measurement time.

BELCONTROL

Quick BET	igen	Multi-point BET surface area in less than 20 minutes
Helium-Free Measurement	igen	AFSM™2 enables He-free measurements with unmatched accuracy
Adsorption Kinetics	opcionális	Rate of adsorption measurements for diffusion analysis

OPERATION SOFTWARE



When developing the software, we focused on simplifying the operations. The step-by-step instructions allow users of all levels of experience to operate the instrument safely and effectively. Simple Mode needs sample labeling, pretreatment conditions (skippable if externally done) and measurement points; Professional Mode provides detailed configuration of dosing, equilibrium and leak-check options.

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GAS ADSORPTION MEASUREMENTS

EVALUATION SOFTWARE BELMASTER


A sorption isotherm is the relationship between the adsorbed amount and the equilibrium pressure of an adsorptive gas at constant temperature. Gas sorption isotherms provide information about specific surface area, pore size distribution and pore volume.

The specific surface area (SSA) refers to the accessible surface area of the sample and is of great importance in adsorption and catalytic reactions. It can be calculated according to ISO 9277 using the BET method or the Langmuir method. Select the pressure range (multi-point BET) or a measurement point (single-point BET) to calculate the surface area automatically. BELMASTER also supports ISO 9277 Annex C for microporous materials.

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C KONSTANS

AFSM™

The following measurements are conducted using our patented Advanced Free Space Measurement (AFSM™). A reference cell tracks the changes in free-space volume for high reproducibility. A certified reference material CRM-170 (nonporous alpha alumina; certified nitrogen BET SSA: $1.05 \pm 0.05 \text{ m}^2/\text{g}$) checks accuracy at low total surface areas. Although the total surface area is reduced from approx. 2 m^2 to below 0.4 m^2 , the BET SSA remains almost identical.

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

Theory	Interaction Surface and Gas	Adsorbate	Applicable Pore Size Range
BJH, CI, DH, INNES	Kelvin equation (Surface tension and contact angle)	Bulk liquid density	> 2 nm; meso and macropore
HK, SF, CY	Lennard-Jones potential (Interaction and repulsion force)	Bulk liquid density	0.4–2 nm; micropore
NLDFT, GCMC	Statistical thermodynamic model		0.35–500 nm; whole pore range

Microtrac provides evaluation methods for N₂ (77.4 K), Ar (87.3 K) and CO₂ (298 K), using NLDFT / GCMC kernels of slit, cylinder and cage pore models with carbon and metal oxide surface atoms. BELMASTER makes it easy to compare experimental and simulated isotherms, the latter serving as a basis for pore size distribution calculations.

In the following part an example for pore size distribution calculation by BJH method is given:

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
Microtrac's gas adsorption analyzers are used in a variety of fields. These include catalysts, batteries, fibers, polymer materials, zeolite, fuel cells, chemicals, pigments, cosmetics, MOF / PCP, magnetic powders, separation membranes, filters, toners, cement, ceramics, semi-conductors, and many more.

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STANDARD CONSUMABLE GOODS

The standard consumables comprise sample cells, filler rods, filters, O-rings, caps and weighing platforms. Various sizes of sample cells, quick seals, and other optional consumables are also available.

WATER BATH

Easily swap the Dewar with a water bath for measurement temperatures ranging from -10°C to 80°C. A refrigerated / heated circulator is required.


HEATER

In-situ pretreatment of the sample from 50°C to 550°C, directly at the measurement port without transferring the sample cell before analysis. Particularly useful for sensitive samples.

GAS SELECTOR

The gas selector enables the connection of up to 4 adsorptives simultaneously (1 × Helium + 4 × adsorptive).

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microtrac.com/hu/contact-services/quote-request-form/?product=602

MICROTRAC
a **VERDER** company

microtrac.com/hu/termek/gas-adsorption-measurement/surface-area-pore-size-distribution/belsorp-mini-x/

© Microtrac – microtrac.com / Content may be subject to modifications or corrections